**Comfort Check Pacifier Thermometer
with BONUS Pacifier Medicine**

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⚠ WARNING

BEFORE USING THESE PACIFIER PRODUCTS, PLEASE READ CAREFULLY.

- **ALWAYS** supervise your child when using thermometer or medicine dispenser. **NEVER** allow your child to use them alone or without close parental supervision. Keep all medications out of child's reach.
- **NEVER** allow children to play with these products. They are not toys. The thermometer is a sensitive electronic device and may be damaged if not used properly under parental supervision.
- **DO NOT** tie pacifier around a child's neck as it presents a strangulation danger.
- **NOT** intended for extended use as a pacifier.
- Avoid sudden impact and dust to thermometer.
- **FOR HOME USE ONLY.**

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To Operate Pacifier Thermometer

1. Press and release the small button. A full LCD screen will display for 2 seconds. After full display, last temperature measured will appear for 2 seconds. The display panel will flash to indicate that it is ready for insertion into the child's mouth to measure temperature. In the display will appear 188.8° F for about two seconds, and then "L" with the °F symbol flashing.
2. Allow the child to suck on the Pacifier Thermometer for a minimum of 90 seconds. The thermometer will beep for 8 to 10 seconds when it reaches maximum temperature and the °F will stop flashing.
3. Remove the Pacifier Thermometer and read the temperature indicated on display panel. **NOTE:** It is not uncommon for infants to have a normal body temperature slightly higher than the prescribed norm for adults.
4. To turn the thermometer off, press the button switch again. In order to conserve battery life we recommend turning the thermometer off after taking a reading. However the thermometer will automatically shut off after about 10 minutes.
5. If the measurement is interrupted for any reason, repeat the above procedure.
6. H°F will be displayed when the measured temperature is above 109.9° F. Turn off, wait one minute, then turn on and try again. If H°F displays again, discard thermometer.
7. To recall the last temperature, press and hold the power button for at least 3 seconds. The thermometer will run a built-in test. The last temperature will be displayed with the letter "M" under the degree sign. Release the power button and the thermometer will return to the ready-to-measure mode.

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CARE AND CLEANING:

- Examine the pacifier bulb before each use for cracks, tears, tackiness or other signs of wear. Discard the worn unit immediately.
- Do not attempt to pull or twist the nipple apart from the body or detach button head as this will damage the thermometer's sensitive temperature probe and electronics.
- The thermometer may be sterilized by immersing the silicone rubber nipple **ONLY** in sterilizing solution such as alcohol. Immerse the nipple approximately 1 inch / 2.5 cm **ONLY** in the solution. The protective clip-on plastic shield supplied may be used as a container for this purpose.
- The plastic parts of the thermometer may be wiped clean with a soft cloth.
- Use this thermometer for oral temperature measurements only.

CAUTION: The pacifier thermometer should never be completely immersed in sterilizing solution, or sterilized in boiling water, as damage will occur to the electronic components. DO NOT place in microwave.

SPECIFICATIONS:

1. The measuring range is 90.0° F to 109.9° F.
2. The accuracy of the Safety 1st Pacifier Thermometer is $\pm 0.2^\circ$ F, from 98.0° F to 102.0° F.
3. The life of the 1.5V drycell battery sealed in the thermometer is approximately 18 months under normal use conditions. When the display shows "BAT", discontinue use. **NOTE:** The battery is not replaceable. Do not attempt to remove the battery.

BONUS Pacifier Medicine Dispenser

1. Remove backing and fill to prescribed dosage.
2. Snap backing onto pacifier.

Care and Cleaning

Examine the pacifier bulb before each use for cracks, tears, tackiness or other signs of wear. Discard the worn unit immediately. Wash with warm, soapy water and rinse thoroughly before and after each use.

One Year Limited Warranty for Thermometer

Dorel Juvenile Group, Inc. warrants to the original purchaser that this product (Comfort Check Pacifier Thermometer), is free from material and workmanship defects when used under normal conditions for a period of one (1) year from the date of purchase. Should the product contain defects in material or workmanship Dorel Juvenile Group, Inc., will repair or replace the product, at our option, free of charge. Purchaser will be responsible for all costs associated with packaging and shipping the product to Dorel Juvenile Group Consumer Relations Department at the address noted on the front page and all other freight or insurance costs associated with the return. Dorel Juvenile Group will bear the cost of shipping the repaired or replaced product to the purchaser. Product should be returned in its original package accompanied by a proof of purchase, either a sales receipt or other proof that the product is within the warranty period. This warranty is void if the owner repairs or modifies the product or the product has been damaged as a result of misuse.

This warranty excludes any liability other than that expressly stated previously, including but not limited to, any incidental or inconsequential damages.

SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS THAT MAY VARY FROM STATE TO STATE.

©2013 Dorel Juvenile Group. All Rights Reserved. Todos derechos reservados.
www.djgusa.com (800) 544-1108 www.safety1st.com
Made in CHINA. Hecho en CHINA.
Styles and colors may vary. Los estilos y los colores pueden variar.
Distributed by (distribuido por) Dorel Juvenile Group, Inc.,
2525 State St., Columbus, IN 47201-7494
Dorel Distribution Canada, 873 Hodge, St. Laurent, QC H4N 2B1
02/26/13 4358-5511B

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